



ExpressConverter 1.9 Release Note

Express Converter is a file transformation tool that reads microarray data files in a variety of file formats and generates TIGR Array Viewer file (.tav) or TIGR MultiExperiment Viewer file (.mev) as output so that the microarray data can be uploaded to the database with MADAM or analyzed with MIDAS and MeV.

Change and Issue Solved

The ScanArray tab has been modified to be able to convert the .csv file, which is one of the outputs of ScanArray Express 2.0, into mev file format. At the same time, the mapping and sorting algorithms of the converting have been enhanced so that it is unnecessary to have the original Genepix and ImaGene files as the data input.